

APA222Ca02 100µg

Active Interferon Beta (IFNβ)

Organism Species: *Canis familiaris*; Canine (Dog)

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression.

Host: *E. coli*

Residues: Ala21~Asn186

Tags: N-terminal His-tag

Purity: >90%

Traits: Freeze-dried powder

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Buffer Formulation: PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

Original Concentration: 200µg/mL

Applications: Cell culture; Activity Assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 6.2

Predicted Molecular Mass: 23.8kDa

Accurate Molecular Mass: 24kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```
          AMSNDLLRSQ LSSSSLECQE LLLQLNGTTE
YCLKDRINFE IP EEIEKSRQ FQKEDIILIT HEMFQKIFDI FRNISRGTGW
NETTVENLLV KLHWQKEHLE IILEDVKEKE NFTWDNRLL HLKYYLRIV
QYLKAKEYSI CAWTIVQAEI CRNFFFLNIL TDYLQN
```

[ACTIVITY]

Interferon Beta (IFN β) is belongs to type I interferons (IFNs) family which a large subgroup of interferon proteins that help regulate the activity of the immune system. The IFN β proteins are produced in large quantities by fibroblasts. They have antiviral activity that is involved mainly in innate immune response. Two types of IFN β have been described, IFN β 1 (IFNB1) and IFN β 3 (IFNB3). IFN β 1 is used as a treatment for multiple sclerosis as it reduces the relapse rate. To test the effect of IFN β on cell apoptosis, A549 cells were seeded into triplicate wells of 96-well plates at a density of 4,000 cells/well and allowed to attach overnight, then the medium was replaced with various concentrations of recombinant dog IFN β diluted with 5% serum standard DMEM. After incubated for 72h, cells were observed by inverted microscope and cell proliferation was measured by Cell Counting Kit-8 (CCK-8). Briefly, 10 μ l of CCK-8 solution was added to each well of the plate, then the absorbance at 450 nm was measured using a microplate reader after incubating the plate for 1-4 hours at 37 °C . Apoptosis of A549 cells after

incubation with IFN β for 72h observed by inverted microscope was shown in Figure 1. Cell viability was assessed by CCK-8 (Cell Counting Kit-8) assay after incubation with recombinant dog IFN β for 72h. The result was shown in Figure 2. It was obvious that IFN β significantly decreased cell viability of A549 cells. The EC50 of recombinant dog IFN β is about 3 μ g/ml.




Figure 1. Cell apoptosis of A549 cells after stimulated with recombinant dog IFN β .
(A) A549 cells cultured in DMEM, stimulated with 5 ug/ml IFN β for 72h;
(B) Unstimulated A549 cells cultured in DMEM for 72h.




Figure 2. Cell apoptosis of A549 cells after stimulated with recombinant dog IFN β .

[IDENTIFICATION]




Figure 3. Gene Sequencing (extract)




Figure 4. SDS-PAGE

Sample: Active recombinant IFN β , Dog

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.