



rP92281Hu01 100µg Adenosine Kinase (ADK) **Organism: Homo sapiens (Human)** Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

7th Edition (Revised in May, 2013)


[PROPERTIES]

Residues: Lys94~Pro339 (Accession # P55263), with N-terminal His-Tag. Host: E. coli **Purity:** >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% sucrose, 0.01% sarcosyl. Predicted isoelectric point: 6.6 Predicted Molecular Mass: 29.0kDa Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.



Unique product S uperb quality Client favorite Nicest service 🧕 ISO9001:2008; 🧟 ISO13485:2003; $\zeta \in$





[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[<u>SEQUENCES</u>]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSG SEF-KFGEILK RKAAEAHVDA HYYEQNEQPT GTCAACITGD NRSLIANLAA ANCYKKEKHL DLEKNWMLVE KARVCYIAGF FLTVSPESVL KVAHHASENN RIFTLNLSAP FISQFYKESL MKVMPYVDIL FGNETEAATF AREQGFETKD IKEIAKKTQA LPKMNSKRQR IVIFTQGRDD TIMATESEVT AFAVLDQDQK EIIDTNGAGD AFVGGFLSQL VSDKPLTECI RAGHYAASII IRRTGCTFP

Unique product Superb quality Client favorite Nicest service @ ISO9001:2008; @ ISO13485:2003; C \in