

P90014Rb02
Bone Morphogenetic Protein 4 (BMP4)
Organism: *Oryctolagus cuniculus* (Rabbit)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Protein Names: Bone Morphogenetic Protein 4

Synonyms: BMP4

Species: Rabbit

Size: 100 μ g

Source: *Escherichia coli*-derived

Subcellular Location: Secreted, extracellular space,
extracellular matrix.

[PROPERTIES]

Residues: Thr211~Lys392 (Accession # O46576),
with N-terminal His-Tag.

Grade & Purity: >95%, 24kDa as determined by
SDS-PAGE reducing conditions.

Formulation: Supplied as lyophilized form in PBS, pH
7.4, containing 5% sucrose, 0.01% sarcosyl.

Endotoxin Level: <1.0 EU per 1 μ g (determined by
the LAL method).


Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 22.5kDa

Predicted isoelectric point: 9.6

Rabbit BMP4



[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF- TRWETFDVSP AVLRWTREKQ PNHGLAVEVT HFHHTRTHQG
QHVRLSRSLL QGSGDWAQFR PLLVTFGHDG RGHALTRRRR AKRSLKHPQ
RARKKNKNCR RHALYVDFSD VGWNDWIVAP PGYQAFYCHG DCPFPLADHF
NSTNHAIVQT LVNSVNSSIP KACCVPTELS AISMLYLDEY DK