

P94212Mu01
Elastase 3B (ELA3B)
Organism: Mus musculus (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

5th Edition (Revised in January, 2013)

[DESCRIPTION]

Protein Names: Elastase 3B

Synonyms: ELA3B, Cela3b, Ela3

Species: Mouse

Size: 100µg

Source: *Escherichia coli*-derived

[PROPERTIES]

Residues: Val28~Ala267 (Accession # Q9CQ52),
with N-terminal His-Tag.

Grade & Purity: >95%, 29kDa as determined by
SDS-PAGE reducing conditions.

Formulation: Supplied as lyophilized form in PBS, pH
7.4, containing 5% sucrose.

Endotoxin Level: <1.0 EU per 1µg (determined by
the LAL method).


Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 27.4kDa

Predicted isoelectric point: 5.6

Mouse ELA3B kDa



15% SDS-PAGE

[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHSGSEF-VVN GEEAVPHSWP WQVSLQYEKD GSFHHTCGGS LITPDWVLTA
GHCISTSRTY QVVLGEHERG VEEGQEQVIP INAGDLFVHP KWNSMCVSCG NDIALVKLSR
SAQLGDAVQL ACLPPAGEIL PNGAPCYISG WGRNSTNGPL PDKLQQALLP VVDYEHCSRW
NWWGLSVKTT MVCAGGDIQS GCNGDSSGGL NCPADNGTWQ VHGVTSEFVSS
LGCNTRLKPT VFTRVSAFID WIEETIA