

MAB180Po23

Monoclonal Antibody to Clusterin (CLU)

Organism Species: *Sus scrofa*; *Porcine (Pig)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Monoclonal antibody preparation

Host: Mouse

Antibody isotype: IgG1 Kappa

Purification: Protein A + Protein G affinity chromatography

Clone number: D24

Traits: Liquid

Concentration: 1mg/ml

UOM: 200µl

Cross Reactivity: N/A

Applications: WB; IHC; ICC; IP.

[IMMUNOGEN]

Immunogen: Recombinant CLU (Ile229~Glu446 (Accession # Q29549)) expressed in *E.coli*

Accession No.: RPB180Po01

[APPLICATIONS]

Western blotting: 0.5-2µg/mL;

Immunohistochemistry: 5-20µg/mL;

Immunocytochemistry: 5-20µg/mL;

Optimal working dilutions must be determined by end user.

[FORMULATION]

Form & Buffer: Supplied as solution form in 0.01M PBS, pH7.4, containing 0.05% Proclin-300, 50% glycerol.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 4°C for frequent use.

Aliquot and store at -20°C for 24 months.


Stability Test: The thermal stability is described by the loss rate. The loss rate was determined

by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[IDENTIFICATION]



Western Blot; Sample: Porcine Serum
 Primary Ab: 2µg/ml Mouse Anti-Porcine
 CLU Antibody Second Ab: 0.2µg/mL
 HRP-Linked Caprine Anti-Mouse IgG
 Polyclonal Antibody (Catalog:
 SAA544Mu19)



DAB staining on IHC-P; Sample:
 Porcine Liver Tissue; Primary Ab:
 20µg/ml Mouse Anti-Porcine CLU
 Antibody Second Ab: 2µg/mL HRP-
 Linked Caprine Anti-Mouse IgG
 Polyclonal Antibody (Catalog:
 SAA544Mu19)

DAB staining on IHC-P;
 Sample: Porcine Pancreas Tissue;
 Primary Ab: 20µg/ml Mouse Anti-
 Porcine CLU Antibody
 Second Ab: 2µg/mL HRP-Linked
 Caprine Anti-Mouse IgG Polyclonal
 Antibody
 (Catalog: SAA544Mu19)

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.