

P84943Hu01
Peptidase Inhibitor 16 (PI16)
Organism: Homo sapiens (Human)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN DIAGNOSTIC OR THERAPEUTIC PROCEDURES

4th Edition (Revised in August, 2012)

[DESCRIPTION]

Protein Names: Peptidase Inhibitor 16

Synonyms: PI16, CRISP9, PSPBP

Species: Human

Size: 100µg

Source: *Escherichia coli*-derived

Subcellular Location: Membrane; Single-pass type I membrane protein.

[PROPERTIES]

Residues: Leu28~Gly440 (Accession # Q6UXB8), with two N-terminal Tags, His-tag and GST-tag.

Grade & Purity: >95%, 71 kDa as determined by SDS-PAGE reducing conditions.

Formulation: Supplied as liquid form in Phosphate buffered saline(PBS), pH 7.4.

Endotoxin Level: <1.0 EU per 1µg (determined by the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 71.4 kDa

Predicted isoelectric point: 5.6

[PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.





[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MRNKKFELGLEFPNLPYYIDGDVCLTQSMIIRYIADKHNLGGCPKERAIEISMLEGAVLDIRYGVSR IAYS KDFETLKVDFLSKLP EML
KMFEDRLCHKTYLNGDHVTHPDFMLYDALDVVLYMDPMCLDAFPKLVCFKKRIEAIQIDKYLKSSKYIAWPLQGWQATFGGGDHPP
KSDGSTSGSGHHHHHSAGLVPRGTAIGMKETA AAKFERQHMDSPDLGTLEVL FQGPLGSEF-LTD EEKRLMVELH
NLYRAQVSPT ASDMLHMRWD EELAAFAKAY ARQCVWGHNK ERGRRGENLF AITDEGMDVP LAMEEWHHER EHYNLSAATC
SPGQMC GHYT QVVWAKTERI GCGSHFCEKL QGVEETNIEL LVCNYEPPGN VKGKRPYQEG TPCSQCPSGY HCKNSLCEPI
GSPEDAQDLP YLVTEAPSFR ATEASDSRKM GTPSSLATGI PAFLVTEVSG SLATKALPAV ETQAPTSLAT KDPPSMATEA
PPCVTTEVPS ILAAHSLPSL DEEPTVFPKS THVPIPKSAD KVTDKTKVPS RSPENSLDPK MSLTGARELL PHAQEEAEAE
AELPPSSEVL ASVFPAQDKP GELQATLDHT GHTSSKSLPN FPNTSATANA TGGRALALQS SLPGAEGPDK PSVVSGLN SG

