

RPB061Mu01 100µg
Recombinant Protein Disulfide Isomerase (PDI)
Organism Species: *Mus musculus* (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[PROPERTIES]

Residues: Asp212~Arg463 (Accession # P09103),
with two N-terminal Tags, His-tag and GST-tag.

Host: *E. coli*

Subcellular Location: Endoplasmic reticulum
lumen. Melanosome. Cell membrane; Peripheral
membrane protein.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg
(determined by the LAL method).


Formulation: Supplied as lyophilized form in 20mM Tris,
150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT,
0.01% sarcosyl, 5% trehalose, and preservative.

Predicted isoelectric point: 5.7

Predicted Molecular Mass: 60.8kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in ddH₂O.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MSPILGYWKI KGLVQPTRL L LEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSR IA YSKDFETLKV
DFLSKLP EML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFP KLVCFK
KRIEAI PQID KYLKSSKYIA WPLQG WQATF GGGDHPPKSD GSTSGSGHHH HHSAGLVPR
GSTAIGMKET AAKFERQHM DSPDLGTLEV LFQGPLGS-DEGRN NFE G EITKEKLLDF
IKHNQLPLVI EFTEQTAPKI FGGEIKTHIL LFLPKSVSDY DGKLSSFKRA AEGFKGKILF
IFIDSDHTDN QRILEFFGLK KEECPAVRLI TLEEEMTKYK PESDELTAEK ITEFCHRFL E
GKIKPHLMSQ EVPEDWDKQP VKVLVGANFE EVAFDEKKNV FVEFYAPWCG HCKQLAPIWD
KLGETYKDHE NIIIAKMDST ANEVEAVKVH SFPTLKFFPA SADRTVIDYN GER

[REFERENCES]

1. Gong Q.-H., *et al.* (1988) *Nucleic Acids Res.* 16:1203-1203.
2. Mazzarella R.A., *et al.* (1990) *J. Biol. Chem.* 265:1094-1101.
3. Church D.M., *et al.* (2009) *PLoS Biol.* 7:E1000112-E1000112.
4. Park J., *et al.* (2013) *Mol. Cell* 50:919-930.