

**RPK210Hu01 100µg**

**Recombinant A Disintegrin And Metalloproteinase**

**With Thrombospondin 10 (ADAMTS10)**

**Organism Species: Homo sapiens (Human)**

***Instruction manual***

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

## **[ PROPERTIES ]**

**Residues:** Arg239~Pro602

**Tags:** Two N-terminal Tags, His-tag and T7-tag

**Accession:** Q9H324

**Host:** *E. coli*

**Subcellular Location:** Secreted, extracellular space, extracellular matrix.

**Purity:** >95%

**Endotoxin Level:** <1.0EU per 1µg  
(determined by the LAL method).


**Formulation:** Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl.

**Predicted isoelectric point:** 7.8

**Predicted Molecular Mass:** 43.7kDa

**Applications:** SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



## **[ USAGE ]**

Reconstitute in sterile PBS, pH7.2-pH7.4.

## **[ STORAGE AND STABILITY ]**

**Storage: Avoid repeated freeze/thaw cycles.**

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## **[ SEQUENCES ]**

The sequence of the target protein is listed below.

RY VETLVVADKM MVAYHGRRDV EQYVLAIMNI VAKLFQDSSL GSTVNILVTR LILLTEDQPT  
LEITHHAGKS LDSFCKWQKS IVNHSGHGNA IPENGVANHD TAVLITRYDI CIYKNKPCGT  
LGLAPVGGMC ERERSCSVNE DIGLATAFTI AHEIGHTFGM NHDGVGNSCG ARGQDPAKLM  
AAHITMKTNP FVWSSCSRDY ITSFLDSGLG LCLNNRPPRQ DFVYPTVAPG QAYDADEQCR  
FQHGVKSRQC KYGEVCSSELW CLSKSNRCIT NSIPAAEGL CQHTIDKGW CYKRVCVPFG  
SRPEGVDGAW GPWTPWGDCE RTCGGGVSSS SRHCDSRPRT IGGKYCLGER  
RRHRSCNTDD CP