

RPA259Mu01 10 μ g
Recombinant Annexin A5 (ANXA5)
Organism Species: *Mus musculus* (Mouse)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[**PROPERTIES**]

Source: Prokaryotic expression.

Host: *E. coli*

Residues: Met1~Asp319

Tags: N-terminal His-Tag

Tissue Specificity: Liver, Skeletal Muscle, Breast, Brain.

Purity: >98%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% sarcosyl, 5%Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.1

Predicted Molecular Mass: 37.0kDa

Accurate Molecular Mass: 37kDa as determined by SDS-PAGE reducing conditions.

[**USAGE**]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[**STORAGE AND STABILITY**]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MATRGTVTDF PGFDGRADAE VLRKAMKGLG TDEDSILNLL TSRSNAQRQE
IAQEFKTLFG RDLVDDLKSE LTGKFEKLIV AMMKPSRLYD AYELKHALKG
AGTDEKVLTE IIASRTPEEL SAIKQVYEEE YGSNLEDDVV GDTSGYYQRM
LVVLLQANRD PDTAIDDAQV ELDAQALFQA GELKWGTDEE KFITIFGTRS
VSHLRRVFDK YMTISGFQIE ETIDRETSGN LEQLLLAVVK SIRSIPAYLA
ETLYYAMKGA GTDDHTLIRV VVSRSEIDLF NIRKEFRKNF ATSLYSMIKG
DTSGDYKKAL LLLCGGEDD

[IDENTIFICATION]




Figure 1. SDS-PAGE