

RPA798Ra01 200µg
Recombinant Hypoxia Inducible Factor 1 Alpha (HIF1a)
Organism Species: Rattus norvegicus (Rat)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)



### [PROPERTIES]

**Source:** Prokaryotic expression.

Host: E. coli

Residues: Pro216~Val336 Tags: N-terminal His-Tag

Tissue Specificity: Liver, Kidney, Brain.

Subcellular Location: Cytoplasm. Nucleus. Nucleus speckle.

**Purity: >92%** 

**Traits:** Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 1mM DTT, 5% trehalose, 0.01%

sarcosyl and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 5.6

Predicted Molecular Mass: 15.4kDa

Accurate Molecular Mass: 16kDa as determined by SDS-PAGE reducing conditions.

#### [USAGE]

Reconstitute in PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

## [STORAGE AND STABILITY]

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

### [SEQUENCE]

PMTCL VLICEPIPHP SNIEIPLDSK TFLSRHSLDM KFSYCDERIT ELMGYEPEEL LGRSIYEYYH ALDSDHLTKT HHDMFTKGQV TTGQYRMLAK RGGYVWVETQ ATVIYNTKDS QPQCIV

# [ IDENTIFICATION ]




Figure 1. SDS-PAGE