

RPA175Mu01 100µg
Recombinant Interferon Alpha 4 (IFNα4)
Organism Species: *Mus musculus* (Mouse)
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

12th Edition (Revised in Aug, 2016)

[PROPERTIES]

Source: Prokaryotic expression

Host: *E.coli*

Residues: Asp26~Glu186

Tags: N-terminal His Tag

Subcellular Location: Secreted

Purity: > 95%

Traits: Freeze-dried powder

Buffer formulation: 20mM Tris, 150mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT, 0.01% SKL, 5% Trehalose and Proclin300.

Original Concentration: 200µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 8.3

Predicted Molecular Mass: 19.9kDa

Accurate Molecular Mass: 20kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 20mM Tris, 150mM NaCl (pH8.0) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```

                DLPHT YNLGNKRALT VLEEMRRLPP
LSCLKDRKDF GFPLEKVDNQ QIQKAQAILV LRDLTQQILN LFTSKDLSAT
WNATLLDSFC NDHQQQLNDL KACVMQEPPL TQEDSLLAVR TYFHRITVYL
RKKKHSLCAW EVIRAEVWRA LSSSTNLLAR LSEEKE
    
```

[IDENTIFICATION]

ATG AAT CGA CCT GCT CAC ACTT ATA ACC TCG GGA ACA GAG AGG GCT TGA CAG T CCT GGA AG AAT GAG AAG ACT C C C C C C T T T T O C T G C C T G A A G G A C A G G A A G G A T T T T G G A T T C C C C T T G G A G G T G G A T A C C A A C A A T C C A G A A G G C T C A A G C C A T C C T T G T C T A A G A C





Figure . Gene Sequencing (extract)



[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.