

RPF518Hu01 100µg

Recombinant Kinesin Family, Member 18A (KIF18A)

Organism Species: Homo sapiens (Human)

Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

10th Edition (Revised in Jan, 2014)

[PROPERTIES]

Residues: Met1~Lys453

Tags: Two N-terminal Tags, His-tag and T7-tag

Accession: Q8NI77

Host: *E. coli*

Subcellular Location: Cell projection, ruffle.

Cytoplasm. Nucleus.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg

(determined by the LAL method).


Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl.

Predicted isoelectric point: 8.6

Predicted Molecular Mass: 55.2kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[USAGE]

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The sequence of the target protein is listed below.

MSVTEEDLCH HMKVVVRVVRP ENTKEKAAGF HKVVHVVDKH ILVFDPKQEE VSFFHGKKT
NQNVIKKQNK DLKFVFDVAVF DETSTQSEVF EHTTKPILRS FLNGYNCTVL AYGATGAGKT
HTMLGSADEP GVMYLTMLHL YKCMDEIKKEE KICSTAVSYL EVYNEQIRDL LVNSGPLAVR
EDTQKGVVVH GLTLHQPKSS EEILHLLDNG NKNRTQHPTD MNATSSRSA VFQIYLRQQD
KTASINQNR IAKMSLIDLA GSERASTGA KGTRFVEGTN INRLLALGN VINALADSKR
KNQHYPYRNS KLTRLLKDSL GGNCQTIMIA AVSPSSVFYD DTYNTLYAN RAKDIKSSLK
SNVLNVNHI TQYVKICNEQ KAEILLKKEK LKAYEEQKAF TNENDQAKLM ISNPQEKEIE
RFQEILNCLF QNREEIRQEY LKLEMLLKEN ELK