

RPF396Hu01 10µg
Recombinant Spondin 2 (SPON2)
Organism Species: Homo sapiens (Human)
Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[**PROPERTIES**]

Residues: Ser119~Val331 (Accession # Q9BUD6),
with two N-terminal Tags, His-tag and GST-tag.

Host: *E. coli*

Subcellular Location: Secreted, extracellular
space, extracellular matrix.

Purity: >95%

Endotoxin Level: <1.0EU per 1µg
(determined by the LAL method).


Formulation: Supplied as lyophilized form in 20mM Tris,
500mM NaCl, pH8.0, containing 1mM EDTA, 1mM DTT,
0.01% sarcosyl, 5% trehalose, and preservative.

Predicted isoelectric point: 5.6

Predicted Molecular Mass: 55.2kDa

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)



[**USAGE**]

Reconstitute in ddH₂O.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MSPILGYWKI KGLVQPTRLLEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID
GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV
DFLSKLP EML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK
KRIEAIQID KYLKSSKYIA WPLQGQWQATF GGGDHPPKSD GSTSGSGHHH HHSAGLVPR
GSTAIGMKET AAKFERQHM DSPDLGTLLEV LFQGPLGSEF-SV HEVFSAPAVP
SGTGQTS AEL EVQRRHSLVS FVVRIVPSPD WFBGVDSL DL CDGDRWREQA ALDLYPYDAG
TDSGFTFSSP NFATIPQDTV TEITSSSPSH PANSFYYPRL KALPPIARVT LVRLRQSPRA
FIPPAPVLPS RDNEIVDSAS VPETPLDCEV SLWSSWGLCG GHCGR LGTKS RTRYVRVQPA
NNGSPCPELE EEAECVPDNC V

[REFERENCES]

1. Manda R., *et al.* (1999) *Genomics* 61:5-14.
2. Clark H.F., *et al.* (2003) *Genome Res.* 13:2265-2270.
3. Li Y., *et al.* (2009) *EMBO J.* 28:286-297.
4. Liao C.H., *et al.* (2010) *Endocr. Relat. Cancer* 17:99-111.