

APB406Mu01 100µg
Active Lipopolysaccharide Binding Protein (LBP)
Organism Species: *Mus musculus (Mouse)*
Instruction manual

FOR RESEARCH USE ONLY
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Prokaryotic expression.

Host: *E. coli*

Residues: Leu226~Pro438

Tags: N-terminal His-tag

Purity: >95%

Endotoxin Level: <1.0EU per 1µg (determined by the LAL method).

Buffer Formulation: PBS, pH7.4, containing 0.01% SKL, 5%Trehalose .

Original Concentration: 200µg/mL

Applications: Cell culture; Activity Assays.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 9.1

Predicted Molecular Mass: 27.5kDa

Accurate Molecular Mass: 27kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

```
                LGIDY SLVAAPQAKA QVLDVMFKGE
IFNRNHRSPV ATPTPTMSLP EDSKQMVYFA ISDYAFNIAS RYVHQAGYLN
FSITDDMLPH DSGIRLNTKA FRPFTPQIYK KYPDMKLELL GTVVSAPILN
VSPGNLSLAP QMEIEGFVIL PTSAREPVFR LGVVTNVFAS LTFNNSKVTG
MLHPDKAQVR LIESKVGMFN VNLFQAFLNY YLLNSLYP
```

[ACTIVITY]

Lipopolysaccharide Binding Protein (LBP) is a soluble acute-phase protein that binds to bacterial lipopolysaccharide (or LPS) to elicit immune responses by presenting the LPS to important cell surface pattern recognition receptors called CD14 and TLR4. This protein is part of a family of structurally and functionally related proteins, including BPI, plasma cholesteryl ester transfer protein (CETP), and phospholipid transfer protein (PLTP). It has been reported that LBP can enhance LPS-stimulated IL-8 secretion by THP-1 cells. To test the bioactivity of recombinant mouse LBP, THP-1 cells were seeded into 24-well plate at a density of 1×10^6 cells/mL including 1 ug/mL LPS, and treated with certain concentrations (0.15625 ug/mL-10 ug/mL) of rmLBP for 24h and IL-8 levels in the cell supernatant were determined by ELISA. IL-8 levels in the cell supernatant of THP-1 cells increased significantly after stimulated with recombinant mouse LBP have shown in Figure 1.




Figure1. IL-8 levels in the cell supernatant of THP-1 induced by recombinant mouse LBP

[IDENTIFICATION]




Figure 2. SDS-PAGE

Sample: Active recombinant LBP, Mouse

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.