

EPA909Hu61 10µg

Eukaryotic Ceruloplasmin (CP)

Organism Species: *Homo sapiens (Human)*

Instruction manual

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

13th Edition (Revised in Aug, 2023)

[PROPERTIES]

Source: Eukaryotic expression

Host: 293F cell

Residues: MYRMQLLSICIALSLALVTNS+Lys20~Gly1065

Tags: C-terminal His Tag

Subcellular Location: Secreted

Purity: > 90%

Traits: Freeze-dried powder

Buffer formulation: PBS, pH7.4, containing 5% Trehalose.

Original Concentration: 400µg/mL

Applications: Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

Predicted isoelectric point: 5.8

Predicted Molecular Mass: 123.1kDa

Accurate Molecular Mass: 120kDa as determined by SDS-PAGE reducing conditions.

[USAGE]

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.6 mg/mL. Do not vortex.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCE]

MYRMQLLSICIALSLALVTNSKEKHYYIGIETTWDYASDHGEKKLISVDTEHSNIYLQNGPDRI
GRLYKKALYLQYTDETFRTTIEKPVWLGLGPIIKAETGDKVYVHLKNLASRPYTFHSHGITYYK
EHGAIYPDNTTDFQRADDDKVYPGEQYTYMLLATEEQSPGEGDGNCVTRIYHSHIDAPKDI
ASGLIGPLIICKKDSLDKEKEKHIDREFVVMFSVVDENFSWYLEDNIKTYCSEPEKVDKDNEDF
QESNRMYSVNGYTFGSLPGLSMCAEDRVKWYLFMGNEVDVHAAFFHGQALTNKNYRID
TINLFPATLFDAYMVAQNPGEWMLSCQNLNHLKAGLQAFFQVQECNKSSSKDNIRGKHVR
HYIIAAEEIIWNYAPSGIDIFTKENLTAPGSDSAVFFEQGTTRIGGSYKKLVYREYTDASFTNRK
ERGPEEEHLGILGPVIWAEVGDITRVTFHNKGAYPLSIEPIGVRFNKNNEGTYYSPNYPNQRS
SVPPSASHVAPTETFTYEWTVPKEVGPTNADPVCLAKMYSSAVEPTKDI FTGLIGPMKICKKG
SLHANGRQKDVDKEFYLFPTVFDENESLLEDNIRMFTTAPDQVDKEDEDFQESNKMMSM
NGFMYGNQPGLTMCCKGDSVVWYLF SAGNEADVHGIYFSGNTYLWRGERRDTANLFPQTSL
TLHMWPDTEGTFNVECLTTDHYTGGMKQKYTVNQCRREQSEDSTFYLGERTYYIAAVEVEW
DYSPQREWEKELHHLQEQNVSN AFLDKGEFYIGSKYKKVYRQYTDSTFRVPVERKAE EHL
GILGPQLHADVGDKVKIIFKNMATRPYSIHAGVQTESSTVPTPLPGETLTYVWKIPERSGAG
TEDSACIPWAYYSTVDQVKDLYSGLIGPLIVCRRPYLKVFNPRRKLEFALLFLVFDENESWYLD
DNIKTYS DHPEKVNKDDEEFIESNKMHAINGRMFGNLQGLTMHVGVDEVNWYLMGMGNEI
DLHTVHFHGHSHFQYKHRGVYSSDVFDIFPGTYQTLEMFPRTPGIWLLHCHVTDHIHAGMET
TYTVLQNE DTKSG

[IDENTIFICATION]




Figure . Gene Sequencing (extract)




Figure. SDS-PAGE

[IMPORTANT NOTE]

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.