

FOR RESEARCH USE ONLY

NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

Basic information

Product No. : TSI501Mu19

Organism Species : Kidneys (Mouse) *Same name, Different species.*

Strain : Mouse **Sex :** Male

Age : 6-8 weeks **Weight :** 20g~30g

Level : SPF **Disease :** Normal Adult

Source : Kidney **Traits :** Paraffin block

Applications : Paraffin slides for pathologic research: IHC, IF and HE, Masson and other stainings.

BIOSAFETY of the Mouse Kidney Tissue of Normal Adult

Tissues are negative for HIV-1, HBV, HCV, mycoplasma, bacteria, yeast and fungi.

USAGE of the Mouse Kidney Tissue of Normal Adult

Upon receiving the fresh tissues, transfer the tissues to liquid nitrogen, -80° C freezer or directly process the tissues and experiment immediately. Paraffin embedding tissue often needs to be cut into slices before use. Upon receiving the paraffin section, after grilled with 60° C of 3-5 hours, it can be stored at 4° C in the short term, otherwise it must be stored at -20° C. Upon receiving the frozen section, immediately transfer it to -80° C or -20° C freezer.

STORAGE of the Mouse Kidney Tissue of Normal Adult

Fresh tissues to be used within 1 month must be stored at -80° C. Paraffin embedding tissue and paraffin section to be used within 1 month may be stored at room temperature, otherwise it must be stored at -4° C (≤3 month) or -20° C (≤6 months). Frozen section must be stored at -20° C (≤3 month) or -80° C (≤6 months). Avoid repeated freeze-thaw cycles.

FIGURE of the Mouse Kidney Tissue of Normal Adult




Figure 1




Figure 2

Figure 1 Detection of Mouse Kidney Tissue of Normal Adult by HE Staining (Optical microscope, ×200)

Figure 2 Detection of Mouse Kidney Tissue of Normal Adult by HE Staining (Optical microscope, ×400)