

**RPR142Mu01 100µg**  
**Recombinant Epoxide Hydrolase 4 (EPHX4)**  
**Organism Species: *Mus musculus* (Mouse)**  
***Instruction manual***

FOR RESEARCH USE ONLY  
NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

---

13th Edition (Revised in Aug, 2023)

**[ PROPERTIES ]**

**Source:** Prokaryotic expression

**Host:** *E.coli*

**Residues:** Lys36~Asp359

**Tags:** N-terminal His Tag

**Subcellular Location:** Membrane

**Purity:** > 97%

**Traits:** Freeze-dried powder

**Buffer formulation:** PBS, pH7.4, containing 0.01% SKL, 5% Trehalose.

**Original Concentration:** 200µg/mL

**Applications:** Positive Control; Immunogen; SDS-PAGE; WB.

(May be suitable for use in other assays to be determined by the end user.)

**Predicted isoelectric point:** 7.79

**Predicted Molecular Mass:** 41.3kDa

**Accurate Molecular Mass:** 38kDa as determined by SDS-PAGE reducing conditions.

**[ USAGE ]**

Reconstitute in 10mM PBS (pH7.4) to a concentration of 0.1-1.0 mg/mL. Do not vortex.

**[ STORAGE AND STABILITY ]**

**Storage:** Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.


Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. The loss rate is less than 5% within the expiration date under appropriate storage condition.

**[ SEQUENCE ]**

KLLWS IGRAPAQTFR  
RAARANPPAC LNDPSLGTHC YVRIKDSGLR FHYVAAGERG KPLMLLLHGF  
PEFWYSWRHQ LREFKSEYRV VALDLRGYGE SDAPAHQESY KLDCLIADIK  
DILDSLGYSK CVLIGHDWGG MIAWLIAVCY PEMIMKLIVI NFPHPVSVFTE  
YILRHPAQLF RSSFYFFQI PRFPEFMFSI NDFKALKHLF TSQSTGIGRK  
GRQLTTEDLE AYVYVFSQPG ALSGPINHYR NIFSCLPLKH HMTTPTLLL  
WGEEDAFMEV EMAEVTKIYV KNYFRLTILS EGSHWLQQDQ PDIVNGLIWA  
FLKEETRRD

**[ IDENTIFICATION ]**



**[ IMPORTANT NOTE ]**

The kit is designed for research use only, we will not be responsible for any issue if the kit was used in clinical diagnostic or any other procedures.